


GOODWIN CASE:


**Aggregate Removal
Development and GMP
Manufacturing for a
Conjugated Theranostic
Product**



Aggregate Removal Development

- Product is a complex biologics molecule used for diagnostic imaging in oncology. It is in late-stage clinical trials.
- Process requires several operations:
 - pH adjustments
 - Heating
 - Time-sensitive operations
 - Transition metal removal
- Upon scale-up, with these operations lead to the formation of high molecular weight (HMW) aggregates exceeded levels outside product specification
- Several strategies were investigated to selectively remove undesired HMW aggregates:
 - Chromatographic resins
 - Chromatographic membranes
 - Filter materials
 - Salt fractionation
- Filtration provides excellent aggregate removal while adding minimal change and scalability to the current process compared to other methods
- A selective filtration process was developed upon scale-up to produce desired product in high monomeric purity

Comparison of Aggregate Removal Strategies



Aggregate Filtration Upon Scale Up

Scale	% HMW Aggregates	% Monomer	% HMW Aggregates	% Monomer
	Before Filtration		After Filtration	
Small	1.7	97.2	N/A	N/A
Intermediate	13.2	86.3	3.0	96.9
Large	14.6	84.7	2.0	97.9
Product Specification: HMW Aggregates ≤ 5.0% Monomeric purity ≥ 90%				

Conclusion

- A robust and selective filtration process to remove HMW aggregates was developed to generate a complex biologic product in high monomeric and radiochemical purity
- Filtration was implemented in current manufacturing process, requiring minimal process change and can be scaled for larger cGMP manufacturing campaigns
- Several successfully cGMP batches have been carried out to supply materials for late-stage clinical trials.